

Now Available From Hummingbird Med

Disposable Closed Blood Draw and Flush Set for Central Lines

The **SAMMI SET**.... **S**ingle **A**ccess to any CL catheter for blood **M**anagement and fluid **M**anagement providing **I**mproved line care and **S**tandardization for **E**very blood **T**ransfer procedure.

BTF-HM-3 SAMMI Set (Blood Transfer and Flush Set)



Reducing Infection Risk for Central Lines

- Only 1 access to any Central Line to draw Waste, Sample and Flush. (including pre-flush)
- Never open to atmosphere during use, providing for a fully closed system blood draw.
- No in line components required. Less setup time, less line manipulation, no troubleshooting.
- Also provides reduced infection risk in ICU and Non-ICU areas where Central Line blood draws are done. (ie. Medical Specialty Floors, Emergency Dept, Outpatient Infusion areas)
- Fully disposable after each use, reducing infection risk.
- Standardizes Central Line blood draw technique, reducing risk for user error and infection.

Provides for Blood Conservation and Reduced Transfusion Costs

- Waste blood can be returned to the patient in a closed loop system.
- Return of waste blood has been shown to reduce anemia risk in critical patients.
- May reduce transfusions due to minimal blood loss. (only 0.3mL waste discard if waste returned)
- Blood transfusion average cost is \$1,500.

Order Number **BTF-HM-3** 50 per box