


Hummingbird Med



Introduces an improved method for line / catheter access to reduce infection risk when accessing or changing central lines...




DPD-HM-2

The “**Dual Port Device**” ... with Needle-free, Swabable, Luer access, Hi-flow valves for improved line care*

- Easily Disinfected....Straight through flow.....non Hemolytic
- Microbial barrier effectiveness up to 140 activations....7 days
- Capable of up to 300 activations without loss of valve integrity
- 100% flushing clearance at 3mL flush when used for blood draw
- High backpressure capability after 300 activations (44psi)
- High flow rate capability: @ 5 PSI, 917mL per/min

- Provides a closed access port and a closed line change port when used with a central line

Robertsite® Needlefree Swabable Valve



Order Number: DPD-HM-2 10 per box