



Hummingbird Med Devices, Inc.

Product Description	Product Number	Box Quantity
HUMMI Micro Draw Closed Needleless Blood Transfer Device Sterile, Single package, Non DEHP, Latex Free	ABG-HM-1	Box of 50 ea.
1mL Self-Venting Syringe with 15 units balanced Heparin, No Needle, Not for Injection, Sterile	ABG-3601-VS	Box of 100 ea.
3mL Self-Venting Syringe with NO HEPARIN, No Needle, Not for Injection, Sterile	ABG-3603-VS	Box of 100 ea.
Micro T Split-Septum T-Connector 6" Minibore Extension Set with Split-Septum T-Connector, Rotating Male Luer Lock and Non-Removeable Slide Clamp. Approx. Priming volume 0.25mL, Non DEHP, Latex Free, Sterile	NMT8046	Box of 100 ea.
Dual Port Device with two integral silicone valves, Non-DEHP, Latex Free, Sterile	DPD-HM-2	Box of 10 ea.

Umbilical Catheter with Integral Low Volume Blood Sampling System

UAC Closed Blood Draw Kit Description

- * Single lumen, sterile, radiopaque catheter with radius tip
- * Hummi Micro-Draw Blood Transfer Device
- * Micro T Split-Septum T-Connector
- * Dual Port Device with two closed-port silicone valves

All components sterile, latex free, DEHP free, non-toxic

Catalog Number	French Size & cm	Length	UAC Material	Quantity per box	
UAC-HM-50P	5.0 fr.	32cm	Polyurethane	5 ea.	
UAC-HM-50S	5.0 fr.	32cm	Silicone	5 ea.	
UAC-HM-43P	4.3 fr.	32cm	Polyurethane	5 ea.	
UAC-HM-35P	3.5 fr.	32cm	Polyurethane	5 ea.	
UAC-HM-35S	3.5 fr.	32cm	Silicone	5 ea.	
UAC-HM-25P	2.5 fr.	25cm	Polyurethane	5 ea.	
Dual Lumen Configuration					
UAC-HM-50P2	3.5 fr.	20ga. / 20ga.	32cm	Polyurethane	5 ea.
UAC-HM-35P2	3.5 fr.	21ga. / 24ga.	32cm	Polyurethane	5 ea.